


quantumTM

ADDITIVE MFG CERVICAL PLATE

The **QuantumTM Anterior Cervical Plate** is intended for anterior screw fixation to the cervical spine. It is to be used in skeletally mature patients as an adjunct to fusion of the cervical spine (C2 to T1). The system is indicated for use in the temporary stabilization of the anterior spine during the development of cervical spinal fusion.



The **QuantumTM Anterior Cervical Plate** is an anterior cervical plate system made from titanium conforming to ASTM F3001. The Quantum Anterior Cervical Plate is intended for anterior screw fixation of the plate to the cervical spine. The fixation construct consists of a cervical plate that is attached to the vertebral body of the cervical spine with self-drilling and self-tapping bone screws using an anterior approach. Plates are available in a variety of lengths (08 mm – 110 mm), addressing multiple levels of fixation (one to five).

The plate incorporates graft visualization holes on the longitudinal center line for intraoperative visualization. Bone screws, conforming to ASTM F136 are available in three diameters (3.75 mm, 4.25 mm, and 4.75 mm) and a variety of lengths (10 mm – 20 mm).

ORDERING INFORMATION



Implants - Plates

Part Number	Description
PCP-1-10	1 Level 10mm long
PCP-1-12	1 Level 12mm long
PCP-1-14	1 Level 14mm long
PCP-1-16	1 Level 16mm long
PCP-1-18	1 Level 18mm long
PCP-1-20	1 Level 20mm long
PCP-1-22	1 Level 22mm long
PCP-1-24	1 Level 24mm long
PCP-2-25	2 Level 25mm long
PCP-2-26	2 Level 26mm long
PCP-2-29	2 Level 29mm long
PCP-2-32	2 Level 32mm long
PCP-2-35	2 Level 35mm long
PCP-2-38	2 Level 38mm long
PCP-2-41	2 Level 41mm long
PCP-2-44	2 Level 44mm long
PCP-3-38	3 Level 38mm long
PCP-3-39	3 Level 39mm long
PCP-3-40	3 Level 40mm long
PCP-3-41	3 Level 41mm long
PCP-3-43	3 Level 43mm long
PCP-3-46	3 Level 46mm long
PCP-3-49	3 Level 49mm long
PCP-3-52	3 Level 52mm long
PCP-3-55	3 Level 55mm long
PCP-3-58	3 Level 58mm long
PCP-3-61	3 Level 61mm long
PCP-3-64	3 Level 64mm long
PCP-3-67	3 Level 67mm long

Instrumentation

Part Number	Description
ACP-1001T-AWL	Awl
ACP-1001T-PAWL	Punch awl
ACP-1001T-TAP	Tap 3.75mm dia.
ACP-1001T-DRI	Screw Insertter
PCP-1000T-SINF	Single drill guide fixed
PCP-1000T-SINV	Single drill guide variable
PCP-1000T-DOUF	Double drill guide fixed
ACP-1000T-CDRI	Cannulated Screw Insertter
ACP-1000T-LOCK	Lock engagement tool
ACP-DRL-10	Drill 10mm depth
ACP-DRL-12	Drill 12mm depth
ACP-DRL-14	Drill 14mm depth
ACP-DRL-16	Drill 16mm depth
ACP-DRL-18	Drill 18mm depth
ACP-DRL-20	Drill 20mm depth
ACP-1000T-THRD	Fixation Pin Threaded
ACP-1000T-PUSH	Fixation Pin Push
ACP-1001T-BEND	Plate Bender
ACP-1001T-HOLD	Plate Holder
ACP-1000T-500	N-RTCH, Sport Grip 3.7, A-O Pull ADPT W/Spin Cap

Screws - Fixed

Part Number	Description
PCP-375-10F	3.75mm dia. 10mm long fixed screw
PCP-375-12F	3.75mm dia. 12mm long fixed screw
PCP-375-14F	3.75mm dia. 14mm long fixed screw
PCP-375-16F	3.75mm dia. 16mm long fixed screw
PCP-375-18F	3.75mm dia. 18mm long fixed screw
PCP-375-20F	3.75mm dia. 20mm long fixed screw
PCP-425-10F	4.25mm dia. 10mm long fixed screw
PCP-425-12F	4.25mm dia. 12mm long fixed screw
PCP-425-14F	4.25mm dia. 14mm long fixed screw
PCP-425-16F	4.25mm dia. 16mm long fixed screw
PCP-425-18F	4.25mm dia. 18mm long fixed screw
PCP-425-20F	4.25mm dia. 20mm long fixed screw

Screws - Variable

Part Number	Description
PCP-375-10V	3.75mm dia. 10mm long variable screw
PCP-375-12V	3.75mm dia. 12mm long variable screw
PCP-375-14V	3.75mm dia. 14mm long variable screw
PCP-375-16V	3.75mm dia. 16mm long variable screw
PCP-375-18V	3.75mm dia. 18mm long variable screw
PCP-375-20V	3.75mm dia. 20mm long variable screw
PCP-425-10V	4.25mm dia. 10mm long variable screw
PCP-425-12V	4.25mm dia. 12mm long variable screw
PCP-425-14V	4.25mm dia. 14mm long variable screw
PCP-425-16V	4.25mm dia. 16mm long variable screw
PCP-425-18V	4.25mm dia. 18mm long variable screw
PCP-425-20V	4.25mm dia. 20mm long variable screw
PCP-475-10V	4.75mm dia. 10mm long variable screw
PCP-475-12V	4.75mm dia. 12mm long variable screw
PCP-475-14V	4.75mm dia. 14mm long variable screw
PCP-475-16V	4.75mm dia. 16mm long variable screw
PCP-475-18V	4.75mm dia. 18mm long variable screw
PCP-475-20V	4.75mm dia. 20mm long variable screw

Cases

Part Number	Description
PCP-1000C-101	Quantum Base
PCP-1000C-102	Quantum Base Lid
PCP-1000C-103	Quantum Case Insert
PCP-1000C-104	Quantum Plate Caddy
PCP-1000C-105	Quantum Plate Caddy Lid
PCP-1000C-106	Quantum Screw Caddy
PCP-1000C-107	Quantum Screw Caddy Lid
PCP-1000C-108	Quantum 4-5 Level Plate Caddy
PCP-1000C-109	Quantum 4-5 Level Plate Caddy Lid



Nvision Biomedical Technologies

4590 Lockhill Selma
San Antonio, TX 78249

☎ 210-545-3713 ext. 112

✉ orders@nvisionbiomed.com

🌐 www.nvisionbiomed.com



Scan to visit
our website.



Scan to view Quantum in
augmented reality.