

# *javelin*<sup>®</sup>

*TAILOR'S BUNION FIXATION SYSTEM*

The Javelin<sup>®</sup> Tailor's Bunion Fixation System is a single-use bone fixation device intended to be permanently implanted for the correction of Tailor's Bunion.



40mm Length

30mm Length

The Javelin System provides a streamlined, comprehensive solution for Tailor's Bunion corrections. The system consists of an intramedullary titanium alloy plate and screws that provide stable fixation for 5th metatarsal bone fragments following a transverse osteotomy.

The plate's low-profile design enables the user to achieve maximum medial correction of the capital fragment. A targeting jig is provided in the system to allow for easy insertion of bi-cortical screws targeting the plate holes in the intramedullary canal.



# ORDERING INFORMATION

## Plates

Part Number	Description
TB-30-4	30mm long plate
TB-40-4	40mm long plate

## Locking Screws

Part Number	Description
TB-24-06L	2.40mm dia. 6mm long
TB-24-08L	2.40mm dia. 8mm long
TB-24-10L	2.40mm dia. 10mm long
TB-24-12L	2.40mm dia. 12mm long
TB-24-14L	2.40mm dia. 14mm long
TB-24-16L	2.40mm dia. 16mm long
TB-30-06L	3.00mm dia. 6mm long
TB-30-08L	3.00mm dia. 8mm long
TB-30-10L	3.00mm dia. 10mm long
TB-30-12L	3.00mm dia. 12mm long
TB-30-14L	3.00mm dia. 14mm long
TB-30-16L	3.00mm dia. 16mm long

## Non-locking Screws

Part Number	Description
TB-24-06N	2.40mm dia. 6mm long
TB-24-08N	2.40mm dia. 8mm long
TB-24-10N	2.40mm dia. 10mm long
TB-24-12N	2.40mm dia. 12mm long
TB-24-14N	2.40mm dia. 14mm long
TB-24-16N	2.40mm dia. 16mm long
TB-30-06N	3.00mm dia. 6mm long
TB-30-08N	3.00mm dia. 8mm long
TB-30-10N	3.00mm dia. 10mm long
TB-30-12N	3.00mm dia. 12mm long
TB-30-14N	3.00mm dia. 14mm long
TB-30-16N	3.00mm dia. 16mm long

## Instruments

Part Number	Description
TB-1000T-DR24	Drill for 2.4mm screw
TB-1000T-DR30	Drill for 3.0mm screw
TB-1000T-JG30	30mm Targeting Jig Drill Guide
TB-1000T-JG40	40mm Targeting Jig Drill Guide
TB-1000T-JGIN	Targeting Jig Plate Inserter
TB-1000T-BRCH	Broach
TB-1001T-DPG	Depth Gauge
TB-1000T-ISRT	Screw Inserter
TB-1000T-MLT	Slotted Mallet
TB-1000T-PLX	Plate Extractor
TB-1000T-BRX	Broach Extractor
CSRW-1000T-320	Kwire Trocar 1.0mm x 150mm
HAN-1000T-AOSM	AO Handle, Small

