

vertex[®]

COMPRESSION STAPLE SYSTEM

The Vertex[®] Nitinol Staple is made of superelastic nitinol and is designed to provide compression to an osteotomy or bone fracture.

The unique design of the staple combined with the mechanical properties of the nitinol, allow the staple to have a great performance offering a 63.5% greater pullout strength than its competitors*.



- Offers a user friendly and ergonomic design.
- Has a low profile, high strength, and easy insertion.
- The system consists of the staple and instrumentation such as drills and drill guides that are straightforward and ergonomic allowing for a quick and smooth staple insertion.
- The system comes terminally sterilized with disposable instrumentation.
- The bridge design and super-elastic nitinol allow for maximum compression.

*Data on file



ORDERING INFORMATION



Symmetrical Leg Staples


Part Number	Description
STK-08-08-08	8mm bridge staple, 8mm x 8mm legs
STK-10-10-10	10mm bridge staple, 10mm x 10mm legs
STK-12-12-12	12mm bridge staple, 12mm x 12mm legs
STK-14-14-14	14mm bridge staple, 14mm x 14mm legs
STK-16-14-14	16mm bridge staple, 14mm x 14mm legs
STK-16-16-16	16mm bridge staple, 16mm x 16mm legs
STK-18-14-14	18mm bridge staple, 14mm x 14mm legs
STK-18-18-18	18mm bridge staple, 18mm x 18mm legs
STK-20-14-14	20mm bridge staple, 14mm x 14mm legs
STK-20-18-18	20mm bridge staple, 18mm x 18mm legs
STK-20-20-20	20mm bridge staple, 20mm x 20mm legs



Asymmetrical Leg Staples*

Part Number	Description
STK-10-12-14	10mm bridge staple, 12mm x 14mm legs
STK-10-14-16	10mm bridge staple, 14mm x 16mm legs
STK-10-16-18	10mm bridge staple, 16mm x 18mm legs
STK-18-14-16	18mm bridge staple, 14mm x 16mm legs
STK-18-16-18	18mm bridge staple, 16mm x 18mm legs

*Asymmetrical legs are special order, please allow 12 weeks. Please contact us for more information.



Instruments

Part Number	Description
STK-1000C-001	Vertex Reusable Tray Base
STK-1000C-002	Vertex Reusable Tray Lid
STK-1000C-003	Vertex Reusable Tray Caddy Base
STK-1000C-004	Vertex Reusable Tray Caddy Lid
STK-1000T-001	Sizing Guide Body 08-10
STK-1000T-002	Sizing Guide Body 12-14
STK-1000T-003	Sizing Guide Body 16-18
STK-1000T-004	Sizing Guide Body 20-25
STK-1000T-007	1.7mm Drill
STK-1000T-008	1.9mm Drill
STK-1000T-009	2.1mm Drill
STK-1000T-010	2.3mm Drill
STK-1000T-011	2.5mm Drill
STK-1000T-DIS1	Staple Distractor
CSRW-1000T-330	1.30mm x 150mm K-Wire



Nvision Biomedical Technologies

4590 Lockhill Selma
San Antonio, TX 78249

☎ 210-545-3713 ext. 112

✉ orders@nvisionbiomed.com

🌐 www.nvisionbiomed.com



Scan to visit
our website.



Scan to view Vertex
in augmented reality.