

healix™

COMPRESSION SCREW SYSTEM

Nvision
biomedical technologies

- The **Healix™** Cannulated Headed and Headless Compression Screw System for Forefoot, Midfoot and Hindfoot applications.
- The **Healix™** Cannulated Forefoot Compression Screw System includes 2.0, 2.5, 3.0 or 3.5 headed and headless screws.
- The **Healix™** Cannulated Midfoot Compression Screw System includes 3.5, 4.0, 4.5 headed and headless screws.
- The **Healix™** Cannulated Hindfoot Compression Screw System includes 5.5, 6.0, 6.5, 7.0 headed and headless screws.
- The **Healix™** Cannulated Compression Screw System sets are modular, supporting surgeon preferences.
- Size specific screw caddies incorporating color coded drills, drill guides, countersinks, drivers, and k-wires.
- **Healix™** Cannulated Compression Screws are Self-Drilling and Self-Tapping.
- Dual Pitch “Double Start” Threads.
- Variable pitch **Healix™** design between proximal and distal threads allows for superior compression without stripping.
- Hexalobal Stay-Drive™ Screw Head Design – Holds the screw on the driver and allows for single hand insertion.

HEALIX™ CANNULATED SCREWS ARE AVAILABLE IN THE FOLLOWING SIZES

FOREFOOT SCREW SYSTEM

2.0 mm Screw Sizes	2.5 mm Screw Sizes	3.0 mm Screw Sizes	3.5 mm Screw Sizes
10 mm	10 mm	10 mm	16 mm
12 mm	12 mm	12 mm	18 mm
14 mm	14 mm	14 mm	20 mm
16 mm	16 mm	16 mm	22 mm
18 mm	18 mm	18 mm	24 mm
20 mm	20 mm	20 mm	26 mm
22 mm	22 mm	22 mm	28 mm
24 mm	24 mm	24 mm	30 mm
26 mm	26 mm	26 mm	32 mm
28 mm	28 mm	28 mm	34 mm
30 mm	30 mm	30 mm	36 mm
		32 mm	38 mm
		34 mm	40 mm
		36 mm	
		38 mm	
		40 mm	

MIDFOOT SCREW SYSTEM

3.5 mm Screw Sizes	4.0 mm Screw Sizes	4.5 mm Screw Sizes
16 mm	18 mm	20 mm
18 mm	20 mm	22 mm
20 mm	22 mm	24 mm
22 mm	24 mm	26 mm
24 mm	26 mm	28 mm
26 mm	28 mm	30 mm
28 mm	30 mm	32 mm
30 mm	32 mm	34 mm
32 mm	34 mm	36 mm
34 mm	36 mm	38 mm
36 mm	38 mm	40 mm
38 mm	40 mm	42 mm
40 mm	42 mm	44 mm
	44 mm	46 mm
	46 mm	48 mm
	48 mm	50 mm
	50 mm	52 mm
		54 mm
		56 mm
		58 mm
		60 mm

HINDFOOT SCREW SYSTEM

5.5 mm Screw Sizes	6.0 mm Screw Sizes	6.5 mm Screw Sizes	7.0 mm Screw Sizes
24 mm	26 mm	28 mm	30 mm
26 mm	28 mm	30 mm	32 mm
28 mm	30 mm	32 mm	34 mm
30 mm	32 mm	34 mm	36 mm
32 mm	34 mm	36 mm	38 mm
34 mm	36 mm	38 mm	40 mm
36 mm	38 mm	40 mm	42 mm
38 mm	40 mm	42 mm	44 mm
40 mm	42 mm	44 mm	46 mm
42 mm	44 mm	46 mm	48 mm
44 mm	46 mm	48 mm	50 mm
46 mm	48 mm	50 mm	52 mm
48 mm	50 mm	52 mm	54 mm
50 mm	52 mm	54 mm	56 mm
52 mm	54 mm	56 mm	58 mm
54 mm	56 mm	58 mm	60 mm
56 mm	58 mm	60 mm	65 mm
58 mm	60 mm	65 mm	70 mm
60 mm	65 mm	70 mm	75mm
65 mm	70 mm	75 mm	80 mm
70 mm	75 mm	80 mm	85 mm
75 mm	80 mm	85 mm	90 mm
	85 mm	90 mm	95mm
	90 mm	100 mm	100 mm
			105 mm
			110 mm
			115mm
			120 mm
			125mm
			130 mm



ORDERING INFORMATION

Nvision Biomedical Technologies

4590 Lockhill Selma Rd

San Antonio, TX 78249

☎ 210-545-3713 ext. 112

✉ orders@nvisionbiomed.com

www.nvisionbiomed.com