



Anterior Cervical Plate System




Tangis® is an anterior cervical plate system for stabilization after cervical fusion. All devices are manufactured from titanium alloy (Ti-6Al-4V ELI) per ASTM136.

The Tangis® Anterior Cervical Plate System is a versatile spinal implant system that provides a full range of plate and screw sizes to ensure a better anatomical fit with little or no contouring required. The Tangis® Anterior Cervical Plate System is intended for anterior interbody screw fixation of the cervical spine, and features the innovative positive lock anti-migration system, which secures two screws with one twist. Tangis® features all-through-one instrumentation, which provides a common guide for drilling, tapping and screw placement.

SYSTEM CHARACTERISTICS

- » Slim and low profile (<2 mm)
- » 1-4 level fusion
- » 3 diameters of screws (3.75mm, 4.25mm, and 4.75mm)
- » Self-tapping / Self-drilling / Fixed / Variable screws
- » Screw Lengths from 10mm to 20mm in 2 mm increments

PLATE LEVEL	PLATE LENGTH (MM)									
1	20	22	24	26	28	30	32	34		
2	33	36	39	42	45	48	51	54		
3	53	56	59	62	65	68	71	74	77	
4	60	64	68	72	76	80	84	88	92	



MAR-TAN-002

1 Level Plates

PRODUCT NUMBER	DESCRIPTION
ACP-1-20	1 LEVEL, 20MM LONG PLATE
ACP-1-22	1 LEVEL, 22MM LONG PLATE
ACP-1-24	1 LEVEL, 24MM LONG PLATE
ACP-1-26	1 LEVEL, 26MM LONG PLATE
ACP-1-28	1 LEVEL, 28MM LONG PLATE
ACP-1-30	1 LEVEL, 30MM LONG PLATE
ACP-1-32	1 LEVEL, 32MM LONG PLATE
ACP-1-34	1 LEVEL, 34MM LONG PLATE

2 Level Plates

PRODUCT NUMBER	DESCRIPTION
ACP-2-33	2 LEVEL, 33MM LONG PLATE
ACP-2-36	2 LEVEL, 36MM LONG PLATE
ACP-2-39	2 LEVEL, 39MM LONG PLATE
ACP-2-42	2 LEVEL, 42MM LONG PLATE
ACP-2-45	2 LEVEL, 45MM LONG PLATE
ACP-2-48	2 LEVEL, 48MM LONG PLATE
ACP-2-51	2 LEVEL, 51MM LONG PLATE
ACP-2-54	2 LEVEL, 54MM LONG PLATE

3 Level Plates

PRODUCT NUMBER	DESCRIPTION
ACP-3-53	3 LEVEL, 53MM LONG PLATE
ACP-3-56	3 LEVEL, 56MM LONG PLATE
ACP-3-59	3 LEVEL, 59MM LONG PLATE
ACP-3-62	3 LEVEL, 62MM LONG PLATE
ACP-3-65	3 LEVEL, 65MM LONG PLATE
ACP-3-68	3 LEVEL, 68MM LONG PLATE
ACP-3-71	3 LEVEL, 71MM LONG PLATE
ACP-3-74	3 LEVEL, 74MM LONG PLATE
ACP-3-77	3 LEVEL, 77MM LONG PLATE

4 Level Plates

PRODUCT NUMBER	DESCRIPTION
ACP-4-60	4 LEVEL, 60MM LONG PLATE
ACP-4-64	4 LEVEL, 64MM LONG PLATE
ACP-4-68	4 LEVEL, 68MM LONG PLATE
ACP-4-72	4 LEVEL, 72MM LONG PLATE
ACP-4-76	4 LEVEL, 76MM LONG PLATE
ACP-4-80	4 LEVEL, 80MM LONG PLATE
ACP-4-84	4 LEVEL, 84MM LONG PLATE
ACP-4-88	4 LEVEL, 88MM LONG PLATE
ACP-4-92	4 LEVEL, 92MM LONG PLATE

3.75mm Diameter Screws


PRODUCT NUMBER	DESCRIPTION
ACP-375-10F	3.75MM X 10MM FIXED SCREW
ACP-375-12F	3.75MM X 12MM FIXED SCREW
ACP-375-14F	3.75MM X 14MM FIXED SCREW
ACP-375-16F	3.75MM X 16MM FIXED SCREW
ACP-375-18F	3.75MM X 18MM FIXED SCREW
ACP-375-20F	3.75MM X 20MM FIXED SCREW
ACP-375-10V	3.75MM X 10MM VARIABLE SCREW
ACP-375-12V	3.75MM X 12MM VARIABLE SCREW
ACP-375-14V	3.75MM X 14MM VARIABLE SCREW
ACP-375-16V	3.75MM X 16MM VARIABLE SCREW
ACP-375-18V	3.75MM X 18MM VARIABLE SCREW
ACP-375-20V	3.75MM X 20MM VARIABLE SCREW

4.25mm Diameter Screws

PRODUCT NUMBER	DESCRIPTION
ACP-425-10F	4.25MM X 10MM FIXED SCREW
ACP-425-12F	4.25MM X 12MM FIXED SCREW
ACP-425-14F	4.25MM X 14MM FIXED SCREW
ACP-425-16F	4.25MM X 16MM FIXED SCREW
ACP-425-18F	4.25MM X 18MM FIXED SCREW
ACP-425-20F	4.25MM X 20MM FIXED SCREW
ACP-425-10V	4.25MM X 10MM VARIABLE SCREW
ACP-425-12V	4.25MM X 12MM VARIABLE SCREW
ACP-425-14V	4.25MM X 14MM VARIABLE SCREW
ACP-425-16V	4.25MM X 16MM VARIABLE SCREW
ACP-425-18V	4.25MM X 18MM VARIABLE SCREW
ACP-425-20V	4.25MM X 20MM VARIABLE SCREW

4.75mm Diameter Screws

PRODUCT NUMBER	DESCRIPTION
ACP-475-10V	4.75MM X 10MM VARIABLE SCREW
ACP-475-12V	4.75MM X 12MM VARIABLE SCREW
ACP-475-14V	4.75MM X 14MM VARIABLE SCREW
ACP-475-16V	4.75MM X 16MM VARIABLE SCREW
ACP-475-18V	4.75MM X 18MM VARIABLE SCREW
ACP-475-20V	4.75MM X 20MM VARIABLE SCREW



ORDERING INFORMATION:

Nvision Biomedical Technologies

1350 N Loop 1604 E. Suite 103
San Antonio, TX 78232

Office: 210-545-3713 ext. 112

Email: Orders@nvisionbiomed.com

MAR-TAN-002