



trigon[®]


OSTEOTOMY WEDGE SYSTEM



Cotton



Evans



The Trigon[®] HA Stand-Alone Wedge is made of PEEK-Optima[®] HA Enhanced. The PEEK material is compounded with HA, making it homogeneous throughout the implant.

The Trigon[®] HA Stand-Alone Wedge is provided in configurations for both Evans and Cotton osteotomies. These configurations are available in a variety of sizes, allowing for precise and reliable correction that accommodate a large patient population.

The Trigon[®] HA Stand-Alone Wedge does not require any additional fixation, such as plating.

Ergonomic instrumentation allow for accurate wedge insertion and precise screw placement.



Nvision[™]
biomedical technologies

MAR-TRG-004 REV C

ORDERING INFORMATION



Evans Implants

Part Number	Description
EVP-16-18-06	PEEK HA Small Evans Wedge 6mm
EVP-16-18-08	PEEK HA Small Evans Wedge 8mm
EVP-16-18-10	PEEK HA Small Evans Wedge 10mm
EVP-16-18-12	PEEK HA Small Evans Wedge 12mm
EVP-18-20-06	PEEK HA Medium Evans Wedge 6mm
EVP-18-20-08	PEEK HA Medium Evans Wedge 8mm
EVP-18-20-10	PEEK HA Medium Evans Wedge 10mm
EVP-18-20-12	PEEK HA Medium Evans Wedge 12mm
EVP-20-22-06	PEEK HA Large Evans Wedge 6mm
EVP-20-22-08	PEEK HA Large Evans Wedge 8mm
EVP-20-22-10	PEEK HA Large Evans Wedge 10mm
EVP-20-22-12	PEEK HA Large Evans Wedge 12mm

Cotton Implants

Part Number	Description
COP-16-16-06	PEEK HA Small Cotton Wedge 6mm
COP-16-16-07	PEEK HA Small Cotton Wedge 7mm
COP-16-16-08	PEEK HA Small Cotton Wedge 8mm
COP-18-20-06	PEEK HA Large Cotton Wedge 6mm
COP-18-20-07	PEEK HA Large Cotton Wedge 7mm
COP-18-20-08	PEEK HA Large Cotton Wedge 8mm

Locking Screws


Part Number	Description
WSC-25-10	2.50mm dia. 10mm long
WSC-25-12	2.50mm dia. 12mm long
WSC-25-14	2.50mm dia. 14mm long
WSC-25-16	2.50mm dia. 16mm long
WSC-25-18	2.50mm dia. 18mm long
WSC-25-20	2.50mm dia. 20mm long
WSC-25-22	2.50mm dia. 22mm long
WSC-25-24	2.50mm dia. 24mm long
WSC-25-26	2.50mm dia. 26mm long
WSC-25-28	2.50mm dia. 28mm long
WSC-25-30	2.50mm dia. 30mm long

Instruments

Part Number	Description
WD-1001T-EVI	Evans Wedge Insertter
WD-1001T-COI	Cotton Wedge Insertter
WD-1000T-SREW	Trigon Screw Insertter
WD-1000T-TRI	Trigon Trial Insertter
WD-1000T-DRL	Drill for 2.5mm Screw
CSRW-1000T-350	K-Wires (2.4mmx230mm)
CSRW-1000T-800	Pickups
CSRW-1000T-501	Small AO Handle
INST-1000T-SPR	Lamina Spreaders
INST-1000T-HINT	Hintermann Pin Distractor
WD-1000T-GRFB	Trigon Graft Packing Block
WD-1000T-GRFT	Trigon Graft Packing Tamp

Trials

Part Number	Description
WD-161806T-EV	Small Evans 6mm Trial
WD-161808T-EV	Small Evans 8mm Trial
WD-161810T-EV	Small Evans 10mm Trial
WD-161812T-EV	Small Evans 12mm Trial
WD-182006T-EV	Medium Evans 6mm Trial
WD-182008T-EV	Medium Evans 8mm Trial
WD-182010T-EV	Medium Evans 10mm Trial
WD-182012T-EV	Medium Evans 12mm Trial
WD-202206T-EV	Large Evans 6mm Trial
WD-202208T-EV	Large Evans 8mm Trial
WD-202210T-EV	Large Evans 10mm Trial
WD-202212T-EV	Large Evans 12mm Trial
WD-161606T-CO	Small Cotton 6mm Trial
WD-161607T-CO	Small Cotton 7mm Trial
WD-161608T-CO	Small Cotton 8mm Trial
WD-182006T-CO	Large Cotton 6mm Trial
WD-182007T-CO	Large Cotton 7mm Trial
WD-182008T-CO	Large Cotton 8mm Trial



Nvision Biomedical Technologies

4590 Lockhill Selma
San Antonio, TX 78249

☎ 210-545-3713 ext. 112

✉ orders@nvisionbiomed.com

🌐 www.nvisionbiomed.com



Scan to visit
our website.



View Evans in
augmented reality.



View Cotton in
augmented reality.

MAR-TRG-004 REV C
06/08/2022