

*t*rigon[®]


CC WEDGE SYSTEM

The Trigon[®] HA Stand-Alone Wedge is made of PEEK Optima[®] HA Enhanced, where the HA is fully integrated throughout the PEEK.

The Trigon[®] HA Stand-Alone Wedge is provided in a configuration for Calcaneocuboid Arthrodesis. This configuration is available in a variety of sizes, allowing for precise and reliable correction that accommodate a large patient population.

Ergonomic instrumentation allow for accurate wedge insertion and precise screw placement.

- Fusion Grade PEEK is a unique compound of PEEK (80%) and hydroxyapatite (20%) for radiolucent stable fixation with osteoconductive and osteoinductive properties.
- A cut guide with depth markings is included to trim nails as needed.
- Tantalum markers at the distal end of each nail provide a depth indicator on imaging.
- Cannulated for k-wire guidance for enhanced surgical accuracy.
- Nails are impacted across application sites to provide solid fixation.



Ordering Information



Calcaneocuboid Implants

Part Number	Description
CC-2525-08	PEEK HA CC Wedge 8mm
CC-2525-10	PEEK HA CC Wedge 10mm
CC-2525-12	PEEK HA CC Wedge 12mm
CC-2525-14	PEEK HA CC Wedge 14mm

Locking Screws

Part Number	Description
WSC-25-10	2.50mm dia. 10mm long
WSC-25-12	2.50mm dia. 12mm long
WSC-25-14	2.50mm dia. 14mm long
WSC-25-16	2.50mm dia. 16mm long
WSC-25-18	2.50mm dia. 18mm long
WSC-25-20	2.50mm dia. 20mm long
WSC-25-22	2.50mm dia. 22mm long
WSC-25-24	2.50mm dia. 24mm long
WSC-25-26	2.50mm dia. 26mm long
WSC-25-28	2.50mm dia. 28mm long
WSC-25-30	2.50mm dia. 30mm long

Instruments

Part Number	Description
WD-1001T-CCI	CC Wedge Inserter
WD-1000T-SREW	Trigon Screw Inserter
WD-1000T-TRI2	Trigon Trial Inserter2
WD-1000T-DRL	Drill for 2.5mm Screw
CSRW-1000T-352	K-Wires (2.4mmx150mm)
CSRW-1000T-800	Pickups
HAN-1000T-AOSM	Small AO Handle
INST-1000T-HRTD	Heart Distractor
INST-1000T-HINT	Hintermann Pin Distractor
WD-1000T-CGPB	Trigon CC Graft Packing Block

Trials

Part Number	Description
WD-2508T-CC	CC Wedge 8mm Trial
WD-2510T-CC	CC Wedge 10mm Trial
WD-2512T-CC	CC Wedge 12mm Trial
WD-2514T-CC	CC Wedge 14mm Trial




Nvision Biomedical Technologies

12449 Silicon Drive
San Antonio, TX 78249

☎ 210-545-3713 ext. 112

✉ orders@nvisionbiomed.com

🌐 www.nvisionbiomed.com



Scan to visit
our website.