

*t*rigon[®]

LAPIDUS WEDGE SYSTEM

The Trigon[®] HA Wedge Fixation System is made of PEEK-Optima[®] HA Enhanced, where the HA is fully integrated throughout the PEEK.

The Trigon[®] Lapidus Wedge is indicated for a first metatarsal-cuneiform lengthening arthrodesis. The wedge is available in a variety of sizes, allowing for precise and reliable correction that accommodate a large patient population.

The Trigon[®] Lapidus Wedge requires ancillary fixation not provided in the Trigon system. Ergonomic instrumentation allow for accurate wedge insertion and precise screw placement.

- Fusion Grade PEEK is a unique compound of PEEK (80%) and hydroxyapatite (20%) for radiolucent stable fixation with osteoconductive and osteoinductive properties.
- A cut guide with depth markings is included to trim nails as needed.
- Tantalum markers at the distal end of each nail provide a depth indicator on imaging.
- Cannulated for k-wire guidance for enhanced surgical accuracy.
- Nails are impacted across application sites to provide solid fixation.



Ordering Information



Lapidus Implants

Part Number	Description
LDP-2405-XXXX	PEEK HA Small Lapidus Wedge 5mm
LDP-2408-XXXX	PEEK HA Small Lapidus Wedge 8mm
LDP-2410-XXXX	PEEK HA Small Lapidus Wedge 10mm
LDP-2412-XXXX	PEEK HA Small Lapidus Wedge 12mm
LDP-2805-XXXX	PEEK HA Medium Lapidus Wedge 5mm
LDP-2808-XXXX	PEEK HA Medium Lapidus Wedge 8mm
LDP-2810-XXXX	PEEK HA Medium Lapidus Wedge 10mm
LDP-2812-XXXX	PEEK HA Medium Lapidus Wedge 12mm
LDP-3205-XXXX	PEEK HA Large Lapidus Wedge 5mm
LDP-3208-XXXX	PEEK HA Large Lapidus Wedge 8mm
LDP-3210-XXXX	PEEK HA Large Lapidus Wedge 10mm
LDP-3212-XXXX	PEEK HA Large Lapidus Wedge 12mm

Locking Screws

Part Number	Description
WSC-25-10	2.50mm dia. 10mm long
WSC-25-12	2.50mm dia. 12mm long
WSC-25-14	2.50mm dia. 14mm long
WSC-25-16	2.50mm dia. 16mm long
WSC-25-18	2.50mm dia. 18mm long
WSC-25-20	2.50mm dia. 20mm long
WSC-25-22	2.50mm dia. 22mm long
WSC-25-24	2.50mm dia. 24mm long
WSC-25-26	2.50mm dia. 26mm long
WSC-25-28	2.50mm dia. 28mm long
WSC-25-30	2.50mm dia. 30mm long

Instruments

Part Number	Description
WD-1000T-SHFT	Lapidus Wedge Inserter Shaft
WD-1000T-SREW	Trigon Screw Inserter
WD-1000T-TRI	Trigon Trial Inserter
WD-1000T-DRL	Drill for 2.5mm Screw
CSRW-1000T-352	K-Wires (2.4mmx150mm)
CSRW-1000T-330	K-Wires (1.3mmx150mm)
HAN-1000T-AOSM	Small AO Handle
WD-1000T-LGRB	Lapidus Graft Packing Block
WD-1000T-LCGXX	Lapidus Parallel Cut Guide [05-12mm]
WD-1001T-ROTJ	Lapidus Rotational Jig
WD-1001T-WSLV	Lapidus Wedge Screw Sleeve
WD-1001T-INJR	Lapidus Interfrag Screw Targeting Sleeve (R)
WD-1001T-INJL	Lapidus Interfrag Screw Targeting Sleeve (L)
WD-1001T-INRS	Interfrag Screw Targeting Sleeve (R) 24mm
WD-1001T-INLS	Interfrag Screw Targeting Sleeve (L) 24mm
WD-1001T-IMAR	Lapidus IM Angle Reducer
WD-1000T-KWSL	Lapidus Kwire Sleeve
CSRW-1000T-200	Healix Compression Screw Depth Gauge

Trials

Part Number	Description
LDP-2405-XXXX	Small Lapidus Wedge 5mm Trial
LDP-2408-XXXX	Small Lapidus Wedge 8mm Trial
LDP-2410-XXXX	Small Lapidus Wedge 10mm Trial
LDP-2412-XXXX	Small Lapidus Wedge 12mm Trial
LDP-2805-XXXX	Medium Lapidus Wedge 5mm Trial
LDP-2808-XXXX	Medium Lapidus Wedge 8mm Trial
LDP-2810-XXXX	Medium Lapidus Wedge 10mm Trial
LDP-2812-XXXX	Medium Lapidus Wedge 12mm Trial
LDP-3205-XXXX	Large Lapidus Wedge 5mm Trial
LDP-3208-XXXX	Large Lapidus Wedge 8mm Trial
LDP-3210-XXXX	Large Lapidus Wedge 10mm Trial
LDP-3212-XXXX	Large Lapidus Wedge 12mm Trial



Nvision Biomedical Technologies

12449 Silicon Drive
San Antonio, TX 78249

☎ 210-545-3713 ext. 112
✉ orders@nvisionbiomed.com
🌐 www.nvisionbiomed.com



Scan to visit
our website.