

*t*rigon[®]

MIDFOOT WEDGE SYSTEM

The Trigon Midfoot Wedge System is made of Fusion Grade PEEK, where hydroxyapatite (HA) is fully integrated throughout the PEEK.

This system is indicated for a lengthening arthrodesis of the talonavicular or naviculocuneiform joints. Standard and Wide footprints are provided, in small and large sizes, to accommodate a wide patient population, allowing for precise and reliable application.

The Trigon Midfoot Wedge System can be combined with ancillary fixation via the staples provided. A targeting guide is provided for screw fixation if preferred.

- Fusion Grade PEEK is a unique compound of PEEK (80%) and hydroxyapatite (20%) for radiolucent stable fixation with osteoconductive and osteoinductive properties.
- A cut guide with depth markings is included to trim nails as needed.
- Tantalum markers at the distal end of each nail provide a depth indicator on imaging.
- Cannulated for k-wire guidance for enhanced surgical accuracy.
- Nails are impacted across application sites to provide solid fixation.



Ordering Information



Wedges

Part Number	Description
MFP-2017-06C	PEEK Small (20x17mm) Standard Midfoot 6mm Wedge
MFP-2119-06C	PEEK Large (21x19mm) Standard Midfoot 6mm Wedge
MFP-2133-06C	PEEK Small (21x33mm) Wide Midfoot 6mm Wedge
MFP-2235-06C	PEEK Large (22x35mm) Wide Midfoot 6mm Wedge

Staples

Part Number	Description
STK-20-20-20	20mm Bridge Staple, 20mm x 20mm Legs
STK-25-20-20	25mm Bridge Staple, 20mm x 20mm Legs

Trials

Part Number	Description
MF-2017T-06C	PEEK Small (20x17mm) Standard Midfoot 6mm Trial
MF-2119T-06C	PEEK Large (21x19mm) Standard Midfoot 6mm Trial
MF-2133T-06C	PEEK Small (21x33mm) Wide Midfoot 6mm Trial
MF-2235T-06C	PEEK Large (22x35mm) Wide Midfoot 6mm Trial

Instruments

Part Number	Description
WD-1000T-MSHFT	Lesser Metatarsal/Midfoot Wedge Inserter Shaft
WD-1000T-MSSG	Midfoot Wedge Standard Screw Targeting Guide Assembly
WD-1000T-LMKB	Midfoot Wedge Standard Screw Targeting Guide Assembly
WD-1000T-MWSG	Midfoot Wedge Wide Screw Targeting Guide
WD-1000T-LMKB	Midfoot Wedge Wide Screw Targeting Guide
HAN-1000T-AOSM	Handle
WD-1000T-MKSL	Midfoot Wedge K-Wire Sleeve
WD-1000T-MSDG	Midfoot Wedge Staple Drill Guide
WD-1000T-TRI	Trigon Wedge Trial Inserter
KWIRE-13-150S	Healix Compression Screw K-Wire (1.3mm x 150mm)
STK-1000T-030	3.0mm Drill
STK-1000T-011	2.5mm Drill
STK-1000T-DIS1	Staple Distractor
STK-1000T-TMP	Tamp



Nvision Biomedical Technologies

12449 Silicon Drive
San Antonio, TX 78249

☎ 210-545-3713 ext. 112

✉ orders@nvisionbiomed.com

🌐 www.nvisionbiomed.com



Scan to visit
our website.