

interference[®]

SCREW SYSTEM WITH MULTI-DRIVE TECHNOLOGY

The multi-drive Interference Screw system is made of Fusion Grade PEEK, a compound of PEEK and hydroxyapatite, to promote osteoconduction.

The multi-drive technology (inner hexalobe, outer hex, and outer flywheel) insertion features provide increased torque resistance and decrease the risk of stripping.

The Multi-Drive Interference Screw System is intended for soft tissue reattachment, i.e. fixation of ligament and tendon graft tissue and tendon transfers. Specifically:

Foot/Ankle:

Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Mid-foot Reconstruction, Metatarsal Ligament Repair, Flexor Hallucis Longus Transfer for Achilles Tendon Reconstruction, and Flexor Digitorum Longus Transfer for Posterior Tibial Tendon Reconstruction

- Fusion Grade PEEK is a unique compound of PEEK (80%) and hydroxyapatite (20%) for radiolucent stable fixation with osteoconductive and osteoinductive properties.
- A cut guide with depth markings is included to trim nails as needed.
- Tantalum markers at the distal end of each nail provide a depth indicator on imaging.
- Cannulated for k-wire guidance for enhanced surgical accuracy.
- Nails are impacted across application sites to provide solid fixation.



Ordering Information



Tray

Part Number	Description
INTF-1001C-101	Implant Caddy
INTF-1001C-102	Implant Caddy Lid
INTF-1001C-103	Instrument Caddy Part 1
INTF-1001C-104	Instrument Caddy Lid Part 1
INTF-1001C-105	Instrument Caddy Part 2
INTF-1001C-106	Instrument Caddy Lid Part 2

Implants

Part Number	Description
INTF-40-10P	4.0mm Diameter x 10mm long, PEEK
INTF-40-15P	4.0mm Diameter x 15mm long, PEEK
INTF-50-15P	5.0mm Diameter x 15mm long, PEEK
INTF-50-20P	5.0mm Diameter x 20mm long, PEEK
INTF-60-20P	6.0mm Diameter x 20mm long, PEEK
INTF-60-25P	6.0mm Diameter x 25mm long, PEEK

Instruments

Part Number	Description
INTF-1000T-DRL40	4.0 Drill
INTF-1000T-DRL50	5.0 Drill
INTF-1000T-DRL60	6.0 Drill
INTF-1001T-INST40	Telescopic inserter 4.0
INTF-1001T-INST56	Telescopic inserter 5.0/6.0
INTF-1000T-MINS40	Multidrive inserter 4.0
INTF-1000T-MINS56	Multidrive inserter 5.0/6.0
HAN-1000T-AOSJ	Handle
INTF-1000T-PASS	Suture Passer
INTF-1000T-SIZE	Tendon Sizer
CSRW-1000T-260	Drill Guide Handle
CSRW-1000T-284	4.0mm Tissue Protectors
CSRW-1000T-286	5.0mm Tissue Protectors
CSRW-1000T-288	6.0mm Tissue Protectors
TB-1001T-DPG	Depth Gauge
CSRW-1000T-345	K-Wire



Nvision Biomedical Technologies

12449 Silicon Drive
San Antonio, TX 78249

☎ 210-545-3713

✉ orders@nvisionbiomed.com

🌐 www.nvisionbiomed.com



Scan to visit
our website.